

Microcentrifuge Tube

(#60015/60115/61015/62015)

1. Applications

Smooth opening and tight closing are the unique features of SPL's Microcentrifuge Tubes. Microcentrifuge Tubes are produced with high quality and high strength polypropylene to meet the demanded of various applications of molecular biology works.

2. Features

- ❖ External graduation
- ❖ Smooth opening and tight closing
- ❖ Choice of 5 colors: blue, pink, yellow, orange, green (Cat. No. 61015)

3. Functional Test

Cytotoxicity	:	Passed
Non-Pyrogenic	:	Passed
RNase/DNase free	:	Passed
Human DNA free	:	Passed

4. Product Range

Cat. No.	Material	Total Volume	Description	Packaging
60015	PP	1.5 ml	Autoclavable	500/5000
60115	PP	1.5 ml	Sterilized	250/5000
61015	PP	1.5 ml	5 Color	500/5000
62015	PP	1.5 ml	Amber	500/5000

5. ISO 9000 Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2008 Quality Management System (K-QA-Q031478) in Korea.

6. Product Photographs




7. Specifications & Technical Drawing

- Specifications

Dimensions		Material (Resin)	Packaging	
			60015/61015/62015	60115
Diameter:	10.33 mm	Polypropylene (PP)	▪ Pieces/Sleeve: 1	▪ Pieces/Sleeve: 1
Height:	38.18 mm		▪ Pieces/Case: 8	▪ Pieces/Case: 8

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr