

Material Safety Data Sheet

Document no: MSDS-27

Product: **SpermSlow™** **Page 1 / 3**

Date first created: 2006/09/19 **Version 3** **Date revised: 2012/02/20**

1. Identification of Substance/Preparation and of the Company

Product name:	SpermSlow™	Manufacturer:	ORIGIO a/s
Cat. No.:	1094		Knardrupvej 2
Product use:	For slowing down the movement of the sperm to allow for the selection of the most mature, viable spermatozoa for ICSI.		DK-2760 Måløv +45 46 79 02 00

2. Composition / Information on Ingredients

Component	CAS#	Approx. %
Water		>95
Sodium Hyaluronate	9067-32-7	<1
Human Serum Albumin		<1
Sodium Chloride	7647-14-5	<1
Amino Acids		<1
Glucose	492-62-6	<0.1
Sodium Bicarbonate	144-55-8	<0.1
Disodium Hydrogen Phosphate	7558-79-4	<0.1
Potassium Chloride	7447-40-7	<0.1
Sodium Pyruvate	113-24-6	<0.01
Nucleotides		<0.01
Calcium Chloride	10043-52-4	<0.01
Magnesium Sulphate	7487-88-9	<0.01
Vitamins		<0.01
Hydrochloric Acid	7647-01-0	<0.01
Sodium Phosphate	13472-35-0	<0.01
Calcium lactate	5743-47-5	<0.01
Sodium Citrate	6132-04-3	<0.01
Potassium Phosphate	7778-77-0	<0.01
L-Malic acid	97-67-6	<0.01
Sodium Acetate	127-09-3	<0.01
Human Insulin Recombinant	11061-68-0	<0.01
Synthetic cholesterol	57-88-5	<0.01
Gentamicin Sulphate	1405-41-0	<0.01

3. Hazard Identification

The product is classified as non-hazardous according to EU Directive 67/548/EEC.

4. First-Aid Measures

Inhalation:	Not a likely source of exposure.
Ingestion:	Wash out mouth with water. If swallowed consult a physician.
Skin:	Wash with soap and water after each contact.
Eyes:	Flush with copious amounts of water. Assure adequate flushing by separating the eyelids with fingers. If irritation develops, consult a physician.

5. Fire Fighting Measures

Fire:	Not considered to be a fire hazard.
Extinguishing media:	Dry chemical, foam, carbon dioxide or water spray.

ORIGIO a/s

Knardrupvej 2 DK-2760 Måløv Denmark Tel.: +45 4679 0200 Fax: +45 4679 0300
Customer Service Tel.: +45 4679 0202 Fax: +45 4679 0302 Web: www.origio.com Mail: origio@origio.com

Material Safety Data Sheet

Document no: MSDS-27

Product:	SpermSlow™	Page 2 /3
Date first created: 2006/09/19	Version 3	Date revised: 2012/02/20

6. Accidental Release Measures

Absorb on suitable absorbent, such as paper tissue.

7. Handling and Storage

Store in original container at 2-8°C, protected from light. Do not freeze.

8. Exposure Controls/Personal Protection

Avoid contact with skin and eyes. Wash thoroughly after handling.

9. Physiological and Chemical properties

Appearance:	A clear, colourless liquid
Odour:	Odourless
Boiling point:	N/A
Vapour Pressure:	N/A
Specific Gravity:	N/A
Vapour Density:	N/A
Melting Point:	N/A
Solubility H ₂ O:	Soluble
Evaporation Rate:	N/A
Flash Point:	None

10. Stability and Reactivity

Stability:	Stable
Conditions to avoid:	None known
Hazardous Decomposition or by-products:	None known
Hazardous polymerisation:	Not expected to occur

11. Toxicological Information

No information available.

12. Ecological Information

No known ecological damage.

13. Disposal Considerations

Dispose of as medical waste.

14. Transportation Information

Classified as non-hazardous for transport.

Material Safety Data Sheet

Document no: MSDS-27

Product:	SpermSlow™	Page 3 /3
-----------------	-------------------	------------------

Date first created: 2006/09/19	Version 3	Date revised: 2012/02/20
---------------------------------------	------------------	---------------------------------

15. Regulatory Information

The product has been evaluated in accordance with EU Directive 67/548/EEC. The product has been classified as non-hazardous.

16. Other Information

The information contained herein is furnished without warranty of any kind. Users should consider these data only as a supplement to other information gathered by them and must take independent determinations of the suitability and completeness of information from all other sources to assure proper use and disposal of these materials and the safety and health of employees and customers. MediCult a/s shall not be held liable for any damage resulting from handling or from contact with its products.

ORIGIO a/s

Knardrupvej 2 DK-2760 Måløv Denmark Tel.: +45 4679 0200 Fax: +45 4679 0300
Customer Service Tel.: +45 4679 0202 Fax: +45 4679 0302 Web: www.origio.com Mail: origio@origio.com